

HST-MPJ25 Metallographic Specimen Grinding Machine



Picture is for reference only

1.Application:

During the preparation of metallographic specimen, the surface of the cut specimen or the surface of the unprocessed specimen is rough and uneven. In order to improve the quality to the specimen preparation and the working efficiency, the cut section of the specimen may be first rectified on a rectifying machine and then be pre-grounded and polished. MPJ-25 Metallographic Sample Grinding Machine is the fit of the hard work of various technicians and various requirements.

2.Technical Parameters:

Model: HST-MPJ25

Platen speed: 1400r/min

Dimension of Abrasive wheel: 250*30*32mm

Electromotor: 750W, 380V, 50Hz

Dimension: 730*420*360mm

Net weight: 70kg

HST GROUP Co., LTD.

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Tel/Whatsapp: +86-159 1008 1986

www.hssdgroup.com - e-mail: admin@hssdtest.com