

HST-ZDD1 Automatic Dropping Device



Picture just for reference

Application:

In the process of preparation of metallographic sample, pre grinding, polishing and grinding are the indispensable procedures.

The suspension needs to drop in the process of grinding and polishing, so this dropping device is just designed for the automatic dropping of the suspension.

This machine is controlled by single chip microcomputer, and it is output by a precise peristaltic pump. The touch panel displays and controls the input speed.

Specification:

Model: HST-ZDD1

Bottle diameter: Dia. 75mm

Time setting range: 0-99 min

Motor: 24V 9W

Dimension: 100*203*245mm

Weight: 4kg