


HST MPT Semi-automatic Grinding Polishing Head



Application:

The HST-MPT type semi-automatic grinding and polishing head is developed and designed through various market research and extensively listening to the requirements of users. It is an economical product that is very suitable for use in small laboratories with moderate sample preparation.

KS-MPT can be used for single sample preparation. The force acting on each specimen can be controlled individually, so that 1, 2 or 3 specimens can be prepared at a time without the need to use blank or "dummy" specimens.



The KS- MPT type grinding and polishing head can be easily connected with a variety of grinding/polishing machines produced by our company. The operation is simple and fast, and the sample preparation quality is high. An ideal choice for unmanned preparation, saving preparation consumables and improving sample preparation efficiency.

Technical Parameters

Model	MPT
Rotate Speed	50 rpm
Sample Force	0~40N
Sample Capacity	1~3 PIECES
Power Supply	220V, 50Hz
Dimension	183*209*252mm
Weight	13KG

Metallographic consumables

Focus on metallographic equipment and consumables research and development of professional brand

			
Diamond cutting piece	Resin cutting sheet	Metallographic cutting Coolant	CBN cutting piece
			
CM-7 Cold mounting material	Cold/hot mosaic powder	Die for cold Mounting	Mounting material
			
Diamond Grinding Disk	Metallographic sandpaper	MEGULAR	Diamond paste
			
Polishing powder	Polishing lubricant coolant	Metallographic polishing cloth	Polishing liquid

HST GROUP Co., LTD.

NO. 4915, West Jingshi Road, Jinan 250012, PR CHINA.

Tel: +86-159 1008 1986

www.hssdgroup.com - e-mail: admin@hssdtest.com